



Mini Plus

Comfortable Compressor Nebulizer

The Mini Plus is a compressor nebulizer designed for easy portability, especially for people on the move. The compact size, light weight and convenient carrying handle of the Mini Plus facilitate portability and allow users to keep their Mini Plus accessible at all times.

Mini Plus

Comfortable Compressor Nebulizer

Features ▼

- Compact and travel-friendly.
- Cost effective and practical.
- Recessed flow outlet fitting for extra protection.
- Motor-piston drive for reliable, consistent performance.
- Delivers 0.5-5.0 micron particles that can descend deeper into the respiratory tract.



Convenient tube holder design.



Flow outlet fitting, air filter and on/off switch are positioned around the periphery on the front panel and recessed for extra protection.



Delivers aerosolized particles of 0.5-5 microns. Dispenses 3 cc of medication in 10 minutes or less. Reliable output for efficient treatment.



Mouthpiece set : Anti-drool "T" mouthpiece with 6" reservoir tube and 7' supply tube.

Mouthpiece set and 2 masks : Anti-drool "T" mouthpiece with 6" reservoir tube, 7' supply tube and one of each adult & pediatric mask.

Specification

Model	Mini Plus	
Electrical Rating	120 VAC, 60 Hz, 1.2 A	230 VAC, 50 Hz, 0.6 A
(Note: Refer to the rating label on the product)	230 VAC, 60 Hz, 0.6 A	
Maximum Compressor Pressure	≥ 35 psi	≥ 35 psi
Operation Flow rate	Approx. 6.0 lpm	Approx. 5.5 lpm
Nebulization Rate	0.2 - 0.3 ml / min	
Particle Size	0.5 - 5 microns	
Dimensions (W×D×H)	14.0 × 18.8 × 10.0 cm	
	5.5" × 7.4" × 3.9"	
Weight	1.6 kg / 3.5 lb	

Distributed by

Apex Medical Corp.
 No. 9, Min Sheng St., Tu-Cheng,
 New Taipei City, 23679, Taiwan
 T +886 2 2268 5568
 F +886 2 2268 6525
 marketing@apexmedicalcorp.com
 www.apexmedicalcorp.com
 Print-2017/10-500 /All rights reserved
 Specifications may change without notice.

