

APEX VacPlus

portable suction unit



APEX VacPlus

The APEX VAC-series suction pump has been designed to combine function with aesthetic appeal and meet the needs of today's hygiene and infection control requirements.

The VacPlus unit is lightweight, compact and can be used in either AC or DC mode. This unit is available with a rechargeable battery and car adaptor.

- AC/DC Mode
- Remote control option
- Integrated canister holder
- Ergonomic easy to carry handle
- East to adjust pressure regulator

APEX VacPlus

portable suction unit



Electrical rating	AC 100~240V,50/60Hz, DC12V
Power consumption	40W
Maximum pressure	600 mmHg (adjustable)
Maximum flow rate	24 lpm
Net weight	3.3kg
Dimensions	352(L) x 206(W) x 192(H) mm
Sound level	<55 dB
2 year warranty excluding battery	

Specifications are subject to change without notice.

For more information please contact us at:

Air Liquide Healthcare, Stratus House, Unit 1 IDA College Park, Blanchardstown Road North, Blanchardstown, Dublin 15, D15 PEC4.

Tel: +353 (0)18091800 Freephone: 1800 24 02 02 (ROI only) 0800 328 5875 (NI only)

Fax: +353 (0)1829 3966 Email: healthie@airliquide.ie Website: www.airliquidehealthcare.ie