

# APEX Mini-Plus

nebuliser compressor



## APEX Mini-Plus

**The Mini-Plus is a compressor nebuliser designed for easy portability, especially for people on the move.**

Chronic respiratory problems are unpredictable and it is important for people to have their treatment devices close at hand.

The compact size, light weight and convenient carrying handle of the Mini-Plus facilitate portability and allow each person to keep their Mini-Plus accessible at all times.

### Key features

- It's easy to carry
- Updated sleek design makes it easy to pack and store
- Integrated handle makes it convenient for travel and mobility
- Motor-piston drive for reliable, consistent performance

# APEX Mini-Plus

nebuliser compressor



<b>Electrical requirements</b>	230V (50Hz)
<b>Average nebulisation rate</b>	0.25ml/min
<b>Maximum pressure</b>	≥35psi
<b>Flow</b>	6.0 lpm
<b>Particle size</b>	0.5 to 5 microns
<b>Compliance</b>	Electrical safety class II type BF CE 0197
<b>Weight</b>	1.6kg (3.5lb)
<b>Dimensions</b>	140 x 188 x 100mm (5.5 x 7.4 x 3.9in)

*Specifications are subject to change without notice.*

For more information please contact us at:

Air Liquide Healthcare, Stratus House, Unit 1 IDA College Park, Blanchardstown Road North, Blanchardstown, Dublin 15, D15 PEC4.

Tel: +353 (0)18091800 Freephone: 1800 24 02 02 (ROI only) 0800 328 5875 (NI only)

Fax: +353 (0)1829 3966 Email: [healthie@airliquide.ie](mailto:healthie@airliquide.ie) Website: [www.airliquidehealthcare.ie](http://www.airliquidehealthcare.ie)