

# ECONOneb

nebuliser compressor



## ECONOneb

**The ECONOneb is a high-end, versatile nebuliser, suitable for all commonly prescribed medicines for chronic chest conditions. Similar to the TURBOneb, it is continuously rated, therefore suitable for frequent use making it ideal for the hospital environment.**

It is designed for patient appeal and is the smallest and lightest model from the Medix range.

### Key features

- Designed for continuous use making it a reliable nebuliser ideal for clinic and hospital use
- Also suitable for home use
- Nebulises all commonly used respiratory drugs and is suitable for all age groups
- Easy to use controls and carrying handle
- High flow rate which is suitable for nebulising bronchodilators (relievers), corticosteroids (preventers) and antibiotics

# ECONOneb

nebuliser compressor



<b>Electrical requirements</b>	220-240V AC (50Hz) 110-120V (60Hz) AC
<b>Typical dynamic flow and pressure</b>	8 lpm @ 20 psi
<b>Particle size</b>	MMAD 3.0um
<b>Nebulisation time</b>	2.5 to 5.5 minutes for 2.5ml
<b>% Particle &lt; 5um</b>	84.1%
<b>Compliance</b>	Electrical safety class II type BF CE 0120 class IIa MMD 93/42/EEC
<b>Weight</b>	2.7kg (5.95lb)
<b>Dimensions</b>	210 x 185 x 185mm (8.3 x 7.3 x 7.3in)

*Specifications are subject to change without notice.*

For more information please contact us at:

Air Liquide Healthcare, Stratus House, Unit 1 IDA College Park, Blanchardstown Road North, Blanchardstown, Dublin 15, D15 PEC4.

Tel: +353 (0)18091800 Freephone: 1800 24 02 02 (ROI only) 0800 328 5875 (NI only)

Fax: +353 (0)1829 3966 Email: [healthie@airliquide.ie](mailto:healthie@airliquide.ie) Website: [www.airliquidehealthcare.ie](http://www.airliquidehealthcare.ie)