

- Committed to patients
- **Everytime, Everywhere**

InnoSpire Deluxe

compressor nebuliser system



Therapy you can trust

The Philips InnoSpire Deluxe is a reliable and easy to use compressor nebuliser system suitable for nebulising most aerosol medications prescribed for common respiratory diseases.

InnoSpire Deluxe is part of the Philips InnoSpire family of compressor nebuliser systems that provide reliable drug delivery to patients.

Suitable for professional or home use, InnoSpire Deluxe has been designed to easily manage prolonged use in a clinical environment as well as frequent use in the home.

Proven nebuliser technology

Coupled either with the SideStream Reusable or SideStream Plus handset, the InnoSpire Deluxe provides highly efficient aerosol therapy in fast treatment times.

Key features

- Matched with SideStream Reusable or SideStream Plus
- 7LPM driving flow
- Continuously rated for constant use
- Large integral storage
- Handle for easy transport





Reliable aerosol delivery

SideStream Reusable and SideStream Plus have been designed to draw more air in, via an active venturi system, resulting in fast drug delivery and short treatment times.

They also have a unique diamond jet design which helps to minimize wear from frequent use and aids consistent delivery of a high quality aerosol.

SideStream Reusable



SideStream Plus

Equipped with a low resistance valve to deliver fast treatments with less waste.







Specifications	Product	
Performance data	SideStream Plus	SideStream Reusable
MMAD	3.42µm	2.93µm
Fine particle fraction	67%	76.3%
Treatment time (2.5mL charge)	3.68 min	4.12 min
Product specifications	InnoSpire Deluxe	
Electricity supply	230V/50Hz,1Amp	
Compressor max. pressure	360kPa (3.6bar)	
Average flow rate	7LPM@10psig	
Max flow rate	9.8LPM	
Thermal overload protector	Thermal fuse; functioning temperature 150°C	
Size	38 x 17 x 11cm	
Weight	2.2kg	
Sound level	58 ± 3dBA	
Operating temperature/ humidity	10°C to 40°C / 10% to 95% non-condensing	

Ordering information	Part number	
Compressor nebuliser systems - including SideStream Plus kit		
InnoSpire Deluxe – 230V/50Hz – EU plug	1110057	
InnoSpire Deluxe – 230V/50Hz – UK plug	1110064	
Compressor nebuliser systems - including SideStream Reusable kit		
InnoSpire Deluxe – 230V/50Hz – EU plug	1112278	
InnoSpire Deluxe – 230V/50Hz – UK plug	1112279	
InnoSpire Deluxe – 220V/60Hz – UK plug	1110061	
InnoSpire Deluxe – 230V/50Hz – EU plug	1110058	
InnoSpire Deluxe – 230V/50Hz – EU plug	1110065	
Accessories/replacement parts		
SideStream Plus reusable kit	1092001	
SideStream reusable year pack	2010A	
SideStream reusable kit with adult mask	1223A	
SideStream reusable kit with pediatric mask	1224A	
SideStream reusable kit with mouthpiece	1225A	
Air inlet filters (pack of 4)	1113525	

For more information please contact us at:

Air Liquide Healthcare, Stratus House, Unit 1 IDA College Park, Blanchardstown Road North, Blanchardstown, Dublin 15, D15 PEC4. Tel: +353 (0)1809 1800 Freephone: 1800 24 02 02 (ROI only) 0800 328 5875 (NI only)

Fax: +353 (0)1829 3966 Email: healthie@airliquide.ie Website: www.airliquidehealthcare.ie