

- Committed to patients
- : Everytime, Everywhere

TURBOneb

nebuliser compressor



TURBOneb

The TURBOneb is a powerful device with the highest flow rate in our range. It can be used on all commonly prescribed medicines but is specifically designed for nebulising thicker medicines. It is continuously rated, therefore suitable for frequent use, making it ideal for the hospital environment.

Its durable and attractive design incorporates a top carrying handle making it easy to transport. Furthermore, it is designed for ease-of-use and can be used on patients of all age groups.

Key features

- Ideal for Hospital and General Practice use
- Sturdy and reliable
- Remains cool during constant use
- Extremely powerful compressor
- Easy to use
- Easy carrying handle

TURBOneb

nebuliser compressor



Electrical requirements	220-240V AC (50Hz) 110-120V (60Hz) AC
Typical dynamic flow and pressure	10.5 lpm @ 20 psi/1.38 bar
Particle size	MMAD 2.8um ²
Nebulisation time	2 to 5 minutes for 2.5ml
% Particle < 5um	81%
Compliance	Electrical safety class II type BF CE 0120 class IIa MMD 93/42/EEC
Weight	2.7kg (5.95lb)
Dimensions	210 x 185 x 185mm (8.3 x 7.3 x 7.3in)

Specifications are subject to change without notice.

For more information please contact us at:

Air Liquide Healthcare, Stratus House, Unit 1 IDA College Park, Blanchardstown Road North, Blanchardstown, Dublin 15, D15 PEC4. Tel: +353 (0)1 809 1800 Freephone: 1800 24 02 02 (ROI only) 0800 328 5875 (NI only)

Fax: +353 (0)1829 3966 Email: healthie@airliquide.ie Website: www.airliquidehealthcare.ie