

Committed to patients

**Everytime, Everywhere** 

## InogenOne® G3 System

Oxygen. Anytime. Anywhere®



The new Inogen One® G3 delivers the independence of previous Inogen One portable concentrators, but it's smaller, lighter, and makes more oxygen per pound.

It's a portable concentrator light enough to be worn care-free by the active oxygen user. That means there are no tanks to refill and no worry about running out of oxygen.

The Inogen One® G3 is your one solution, your single solution, for oxygen at home or away. The Inogen One® G3 is Oxygen. Anytime. Anywhere®.



## The Inogen One® G3 System (IS-300)

- 45% Smaller
- 30% Lighter
- 27% More Oxygen Per Pound

Compared to the Inogen One G2



Product specifications	
Size	Length: 8.75" Width: 3.0" Height: 7.25"
Weight	4.9 lbs. (includes single battery)
Oxygen flow	Intelligent Delivery Technology Four flow settings 1, 2, 3, 4
Power sensing	AC Power Supply 100-240V, 50-60Hz (auto sensing to allow worldwide use) DC Power Cable: for mobile use in car
Battery	Duration (single battery): 3.5 to 4 hours Recharge approximately 2-4 hours with AC or DC power Duration (double battery): 7 to 8 hours Recharge approximately 4-8 hours with AC or DC power
Noise level	42 Decibels*
Warranty	3 year**
Operation	Simple control functions and easy-to- read LCD display
Use	Designed for 24/7 use, at home or away
* At cotting 2	

- \* At setting 2
- \*\* See Inogen One® Limited Warranty Statement for complete warranty information



Inogen One® AC Power Supply\* Item # BA-301



Inogen One® G3
DC Power Cable\*
Item # BA-306




Inogen One® G3
Carry Bag\*
Item # CA-300



Inogen One® G3
Single Lithium Ion
Battery\*
Up to 4 hours run time
Item #BA-300



Inogen One® G3
Double Lithium Ion
Battery\*\*
Up to 8 hours run time
Item #BA-316



\*\* Accessory sold separately



