

Committed to patients

**Everytime, Everywhere** 

# Invacare® Platinum® 9

reduce cost with Invacare high flow concentrator



The Invacare Platinum 9 offers providers a value added and reliable global solution to treat patients with high flow requirements. Use as a stand alone device or couple together with HomeFill\* for complete patient care.

While optimising oxygen concentration, Platinum 9 minimises costs associated to maintenance and delivery. It allows patients to stay at home comfortably and safely thanks to a series of integrated alarms. Users will appreciate its simplicity of use combined with maximum security.





# Invacare® Platinum® 9

#### High performance and reliability combined with simplicity

Powerful and reliable, Invacare Platinum 9 is also user friendly. With its 9 litre oxygen capacity and 96-87% purity, Platinum 9 is dedicated to patients with high oxygen needs.


Easy to manoeuvre thanks to its integrated handle and casters, Platinum 9 is also easy to use and guarantees extra safety for its user.

#### Value added associated to cost control

Cost linked to inventory is kept to a minimum thanks to the small number of spare parts and specific accessories. With no maintenance for three years, Invacare offers you the best value added concentrator.

### Features and accessories







Only one filter to change when dirty



HomeFill® connection port

## Technical data

Platinum 9













Platinum 9









87-96%





ISO8359 MDD

# SensO<sub>2</sub>-System

#### O<sub>2</sub>-concentration

higher than 85% - green light -Normal

between 73% and 85% - yellow light -Check filter and flow rate

lower than 73% - red light -Call service technician

# Service interval

No scheduled maintenance for 3 years. Replace preventative maintenance items as required.

#### For more information please contact us at: