

Committed to patients

Everytime, Everywhere

DreamStation Auto CPAP



DreamStation helps support long-term patient use while creating **new efficiencies for your business**.

Connectivity options offer more monitoring choices than any leading CPAP brand with integrated Bluetooth[®] on every device, and offering optional Wi-Fi and cellular modems.*

Performance Check simplifies in-home device evaluation and can help reduce the number of normally functioning devices that are returned for service.

EZ-Start's automatic, personalized adjustments to CPAP pressure can help patients gradually acclimate to their prescribed level of therapy.

Daily Progress Feedback displays a simple trend of patients' nightly use to help track progress on the path to adherence.

DreamStation features

- Advanced event detection
- Mask type resistance control
- Mask fit check
- Oximetry capable
- Encore patient management system with optional Adherence Profiler
- DreamMapper patient selfmanagement system compatibility

Humidification

- Fixed, adaptive, heated tube
- Pre-heat humidifier option
- Modular design with one-piece, dishwasher-safe water chamber

Protection programs

- Flex Promise program
- Bi-level Rescue program
- "One-Touch" warranty programs

Auto CPAP features

- C-Flex, C-Flex+, A-Flex
- SmartRamp
- EZ-Start
- CPAP-Check
- Auto-Trial
- Opti-Start
- Digital AutoTrak
- Advanced event detection
- Waveform data
- Standard integrated Bluetooth

*Internal assessment of 2015 competitive CPAP data comparing to ResMed Airsense10/Aircurve10 platform and Fisher & Paykel Icon series platform. The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and the use of such marks by Philips Respironics is under license. Other trademarks and trade names are those of their respective owners.

DreamStation Auto CPAP

Comfort features			Encore data capture and reporting	
Modes	Features	Flex Comfort	Advanced Event Detection	Flow Waveforms
Fixed CPAP CPAP-Check Auto-CPAP	EZ-Start Auto-Trial Opti-Start	C-Flex C-Flex+ A-Flex	Standard	Standard

Device ordering information	Part number
DreamStation Auto CPAP	DSX500S11
DreamStation Auto CPAP with humidifier	DSX500H11
DreamStation Auto CPAP with humidifier and heated tube	DSX500T11

Product specifications

Pressure range	$4 - 20 \text{cm} \text{H}_2\text{O}$
Flex pressure relief	0 to 3
Ramp time	0 to 45 min (5-minute increments)
Starting ramp pressure	4 cm H2O to CPAP Pressure (CPAP); 4cm H ₂ O to Min Pressure (Auto CPAP)
Humidification	Heated humidification: fixed, adaptive, heated tube modes Heated tube (optional)
Data storage capacity (minimum)	SD card: >year On-board: <6 months
Altitude compensation	Automatic
Filters	Reusable pollen and disposable ultra-fine
Device controls	LCD, control dial/push button
Warranty	2 years (US)
Dimensions	15.7 (L) x 19.3 (W) x 8.4 cm (H) (without humidifier) 29.7 (L) x 19.3 (W) x 8.4 cm (H) (with humidifier)
Weight	1.33kg/2.94lbs (without humidifier) 1.98kg/4.37lbs (with humidifier) including power supply
Electrical requirements	100-240 VAC, 50/60 Hz, 2.0-1.0 A

Ordering information	
Humidification	Part number
DreamStation humidifier	DSXH
DreamStation water chamber	1122520
RP-DS Humidifier Dry Box Assembly	1120668
RP-DS Humidifier Flip Lid Seal	1120617
RP-DS Humidifier Dry Box Inlet Seal	1120613
Tubing	Part number
Heated tube (15mm)	HT15
Standard tube (15mm)	PR15
Filters	Part number
Pollen (foam) filter, reusable (1/pack)	1122446
Ultra-fine filter, disposable (1/pack)	1122447
Ultra-fine filter, disposable (2/pack)	1122518
Ultra-fine filter, disposable (6/pack)	1122519
Accessories	Part number
Accessories DreamStation link module	Part number 1120293
DreamStation link module	1120293
DreamStation link module DreamStation SpO ₂ module	1120293 1121694
DreamStation link module DreamStation SpO ₂ module Reuseable finger sensor (clip)	1120293 1121694 936
DreamStation link moduleDreamStation SpO2 moduleReuseable finger sensor (clip)Oximetry finger sensor, Adult, Flex	1120293 1121694 936 953
DreamStation link moduleDreamStation SpO2 moduleReuseable finger sensor (clip)Oximetry finger sensor, Adult, FlexFlexiwrap, sens tape, Adult (25/pack)	1120293 1121694 936 953 954A
DreamStation link module DreamStation SpO2 module Reuseable finger sensor (clip) Oximetry finger sensor, Adult, Flex Flexiwrap, sens tape, Adult (25/pack) Reuseable finger sensor (soft), medium	1120293 1121694 936 953 954A 1062050
DreamStation link module DreamStation SpO ₂ module Reuseable finger sensor (clip) Oximetry finger sensor, Adult, Flex Flexiwrap, sens tape, Adult (25/pack) Reuseable finger sensor (soft), medium	1120293 1121694 936 953 954A 1062050 1062051
DreamStation link module DreamStation SpO ₂ module Reuseable finger sensor (clip) Oximetry finger sensor, Adult, Flex Flexiwrap, sens tape, Adult (25/pack) Reuseable finger sensor (soft), medium Reuseable finger sensor (soft), small	1120293 1121694 936 953 9534 1062050 1062051 Part number
 DreamStation link module DreamStation SpO₂ module Reuseable finger sensor (clip) Oximetry finger sensor, Adult, Flex Flexiwrap, sens tape, Adult (25/pack) Reuseable finger sensor (soft), medium Reuseable finger sensor (soft), small Modems Cellular modem, oximetry capable 	1120293 1121694 936 953 9534 1062050 1062051 Part number 100600C
 DreamStation link module DreamStation SpO₂ module Reuseable finger sensor (clip) Oximetry finger sensor, Adult, Flex Flexiwrap, sens tape, Adult (25/pack) Reuseable finger sensor (soft), medium Reuseable finger sensor (soft), small Cellular modem, oximetry capable Cellular modem, non-oximetry 	1120293 1121694 936 953 954A 1062050 1062051 Part number 100600C
 PreamStation link module DreamStation SpO₂ module Reuseable finger sensor (clip) Oximetry finger sensor, Adult, Flex Flexiwrap, sens tape, Adult (25/pack) Reuseable finger sensor (soft), medium Reuseable finger sensor (soft), small Cellular modem, oximetry capable Cellular modem, non-oximetry 	1120293 1121694 936 953 954A 1062050 1062051 Part number 100600C 100605C 100700W
 PreamStation link module DreamStation SpO₂ module Reuseable finger sensor (clip) Oximetry finger sensor, Adult, Flex Flexiwrap, sens tape, Adult (25/pack) Reuseable finger sensor (soft), medium Reuseable finger sensor (soft), small Cellular modem, oximetry capable Cellular modem, non-oximetry Wi-Fi accessory, non-oximetry 	1120293 1121694 936 953 954A 1062050 1062051 Part number 100600C 100605C 100700W 100705W

For more information please contact us at:

Air Liquide Healthcare, Stratus House, Unit 1 IDA College Park, Blanchardstown Road North, Blanchardstown, Dublin 15, D15 PEC4. Tel: +353 (0)1809 1800 Freephone: 1800 24 02 02 (ROI only) 0800 328 5875 (NI only) Fax: +353 (0)1829 3966 Email: healthie@airliquide.ie Website: www.airliquidehealthcare.ie